

# SciCan STATIM® G4

Cassette Autoclaves



## Sterilization and Dryness at speeds faster than conventional chambered autoclaves



### Fast Cycle Times

STATIM 2000 G4 can sterilize unwrapped loads in under 7 minutes, and wrapped loads in under 11 minutes.

STATIM 5000 G4 can sterilize unwrapped loads in under 9 minutes, and wrapped loads in under 16 minutes.



### Reliable

For over 25 years, STATIM's unmatched speed, reliability, and dependability have made it one of the world's best selling cassette autoclaves.



### Capacity

STATIM G4 autoclaves not only offer fast sterilization times but also enough capacity to meet the needs of a busy practice.

STATIM 2000 G4 capacity allows for up to 4 pouches and the STATIM 5000 G4 allows for up to 10 pouches.



### G4 Technology

SciCan's G4 equipment can connect and send information to your smart devices so you can have cycle data, instructions, and maintenance notifications at your fingertips.



### Gentle

The fresh water feature guarantees fresh water with every cycle. State of the art steam technology is gentle on both solid and hollow instruments.



## STAT/IM 2000 G4 Specifications

Cycle Name	Load	Sterilization Temperature	Sterilization Time	Total Time (mm:ss)*	
				Warm Start	Cold Start
SOLID UNWRAPPED (N)	1.0kg	134°C	3.5 min	6:45	9:15
HOLLOW UNWRAPPED (S)	1.0kg	134°C	3.5 min	8:05	11:45
HOLLOW UNWRAPPED (S)	1.0kg	134°C	18 min	22:35	26:15
HOLLOW WRAPPED (S)**	1.0kg	134°C	3.5 min	10:40	15:35
HOLLOW WRAPPED (S)	1.0kg	134°C	18 min	25:10	30:05
RUBBER/PLASTIC (S)	0.4kg	121°C	15 min	18:40	20:15
RUBBER/PLASTIC (S)	0.4kg	121°C	30 min	33:40	35:15
AIR DRYING ONLY			Not a sterilization cycle		

**Unit size (L x W x H):**  
19.5" x 16.25" x 6" / 49.5 cm x 41.5 cm x 15 cm

**Cassette internal dimensions:**  
11" x 7" x 1.4" / 28 cm x 18 cm x 3.5 cm

**Minimum shelf height for proper leveling:**  
8.5" / 22 cm

**Weight without water:**  
46 lbs / 21 kg

**Power consumption:**  
220-240 V, 50/60Hz, 6A

\* Drying not included  
\*\*Cycle for routine check (STAT/IM PCD test)

## STAT/IM 5000 G4 Specifications

Cycle Name	Load	Sterilization Temperature	Sterilization Time	Total Time (mm:ss)*	
				Warm Start	Cold Start
SOLID UNWRAPPED (N)	1.5kg	134°C	3.5 min	8:45	13:15
HOLLOW UNWRAPPED (S)	1.5kg	134°C	3.5 min	10:50	17:30
HOLLOW UNWRAPPED (S)	1.5kg	134°C	18 min	25:20	32:00
HOLLOW WRAPPED (S)**	1.5kg	134°C	3.5 min	15:30	24:00
HOLLOW WRAPPED (S)	1.5kg	134°C	18 min	30:00	38:30
RUBBER/PLASTIC (S)	0.4kg	121°C	15 min	20:20	22:50
RUBBER/PLASTIC (S)	0.4kg	121°C	30 min	35:20	37:50
AIR DRYING ONLY			Not a sterilization cycle		

**Unit size (L x W x H):**  
23.6" x 16.1" x 7.5" / 60 cm x 41 cm x 19 cm


**Cassette internal dimensions:**  
15" x 7" x 3" / 38 cm x 18 cm x 7.5 cm

**Minimum shelf height for proper leveling:**  
10" / 26 cm

**Weight without water:**  
73 lbs / 33 kg

**Power consumption:**  
220-240 V, 50/60Hz, 6A

\* Drying not included  
\*\*Cycle for routine check (STAT/IM PCD test)



### VistaPure™ Water Filtration

Provides two grades of high-quality water for autoclaves, dental bottles, ultrasonic cleaners and the rinse cycle of instrument washers.

V3000

Manufactured by:  
**SciCan Ltd.**  
1440 Don Mills Rd.  
Toronto, ON, M3B 3P9 / Canada,  
T +1 416-445-1600  
TF +1 800-667-7733  
customerservice@scican.com



### STATMAT™ Protect your STAT/IM

Securely seat cassettes atop your STAT/IM with a safe and hygienic silicone mat.

01-112867S - For STAT/IM 2000  
01-112868S - For STAT/IM 5000

VistaPure™ is manufactured by:  
**Crosstex International, Inc.**  
1733 W. Parkside Lane  
Phoenix, AZ 85027 / USA  
T +1 800-531-2221  
VistaResearchGroup.com



### STAT-DRI™ Plus Drying Agent

Enhance drying in your STAT/IM Cassettes with the STAT-DRI Plus drying agents.

2OZPLUS - Stat-Dri Plus 2 oz with Sprayer  
8OZPLUS - Stat-Dri Plus 8 oz with Sprayer  
8OZPLUS - Stat-Dri Plus 8 oz with Cap (refill)  
32OZPLUS - Stat-Dri Plus 32 oz with Cap (refill)  
T473 - Spray Trigger Only for 32 oz Bottle